

Register Computers

This document provides a step-by-step tutorial on how to Register a Computer as a testing station. Registering computers is a one-time occurrence. A certified coordinator is responsible for completing or overseeing these steps at each site.



1



	Name:	ers\dmontgomery\Do	wnloads\Ca	isasTests(26).e	xe
	Publisher:	Unknown Publisher			
	Type:	Application			
	From:	C:\Users\dmontgomery	ads	\CasasTests(26	j
			Run	Cancel	
i Alw	ays ask before	e opening this file			

 Note! Some browsers present the image displayed at the left. Click Run and wait for the program to load.

STATION MANAGEMENT

First person *initiates* registration

CAS Station Ma	AS eTests nagement Logir
Agency(/Site):	000123
User:	coordinator
Password:	•••••
	Login 🖕

4) Coordinator enters user credentials.

a. Agency ID

- i. Provided by CASAS.
- b. User
 - i. Use e-mail address.
- c. Password
 - i. Provided by Agency.
 - ii. Provided by CASAS *only* for ET Basic and

TE Enhanced users.

5) Click Login

Station Registration				
Site Lab Station Name	Building 1 Building 2			
Registration Status	Not Registered			
Identification Info	MAC: "90-B1-1C-7A-FC-35", Model: "OptiPlex 990", CPU: "Intel(R) Core(TM) i5-2400 CPU @ 3.10GHz", OS: "6.1.7601"			
Register Unregister				

STATION REGISTRATION

- 6) Select Site
 - a. Name, campus, or building location of machine(s).
 - i. Set-up by CASAS from submission of Online Site Agreement.

Register Computers



7) Enter Lab name.

- a. Room location of machine(s).
- b. Established by Agency.

8) Enter Station Name

- a. Standard method for identifying computers.
- b. Established by Agency.
- Note! Station Registration information such as the machine's MAC address is recorded on the hosted server.

9) Click Register



Station Registration





Station Registration						
To confirm registration of this m	achine, another user must log in and click 'Register'.					
Site	Building 1					
Lab	22					
Station Nar Another user must confirm this registration!						
Registration.	ОК					
	MAC: "90-B1-1C-7A-FC-35",					
	"Intel(R) Core(TM) i5-2400 CPU					
Identification Info	@ 3.10GHz", OS: "6.1.7601"					
Regis	ster Unregister					
СЛС	ЛC					
C D e Tests						
·· · · · ·						
Station Management Login						
Agonov//Sito):						
Agency(/Site).	000123					
User:	proctor					
Password:	•••••					
Login						
Statio	n Registration					
To confirm registration of this m	achine, another user must log in and click 'Register'.					
Site	Building 1					
Lab	Assessment Center					
Station Name	1-01					
Registration Status	Pending					
	MAC: "90-B1-1C-7A-FC-35",					
_	Model: "OptiPlex 990", CPU: "Intel(R) Core(TM) i5-2400 CPU					
Identification Info	@ 3.10GHz", OS: "6.1.7601"					
Register Unregister						

Another user must confirm this registration!



Second person validates registration

- **11) Proctor** (or second coordinator) enters user credentials.
 - a. Agency ID
 - i. Provided by CASAS.
 - b. User
 - i. Use e-mail address.
 - c. Password
 - i. Provided by Agency.
 - ii. Provided by CASAS *only* for ET Basic and TE Enhanced users.

12) Click Login

The Station Registration remains "Pending" until the second person validates the registration.

- **13) Proctor** (or second coordinator) clicks Register
- Note! The second person should only click "Register." No registration information should be changed.

Register Computers





Important Note! Confirm that the computer clock (date and time) is synchronized with an Internet time server on all stations. An on-screen warning will display when the program launches if the online system detects a time difference between the station and the Internet time server. For directions on fixing time differences; refer to <u>Clock Synchronization</u> (Appendix C).